


Experiment Setup for ^2H -free samples

Two options for running experiments on samples prepared in protio-solvents: (1) add a capillary with lock-solvent, or (2) run without lock.

(1) Setting lock parameters for use with capillary insert.

Problem: lock parameters are set in the system based on the strength of the lock signal, if a capillary is used there are likely to be problems.



Solution: In Topspin, lock to the solvent in the capillary and run the *loopadj* command which optimises the lock regulation parameters. Save as a new solvent, e.g. c6d6_capillary.

If *loopadj* fails due to insufficient signal intensity, copy the parent solvent to a new solvent and modify the values manually – typically change loop time by -21, loop time by +4 and reduce the loop filter to about 21.

$\text{C}_6\text{D}_6 \rightarrow \text{C}_6\text{D}_6\text{-capillary}$, loop gain 0 \rightarrow -21, loop time 0.2 \rightarrow 0.6, loop filter 300 \rightarrow 20.

DMSO \rightarrow DMSO-capillary, loop gain 5 \rightarrow -17, loop time 0.25 \rightarrow 0.65, loop filter 500 \rightarrow 22.

Experiment Setup for ^2H -free samples

Two options for running experiments on samples prepared in protio-solvents: (1) add a capillary with lock-solvent, or (2) run without lock.

(2) Setting the automation to shim using the biggest ^1H solvent signal.

Problem: the system is setup to lock to a deuterium signal by default, and then shim using that signal, which is not possible if there is no deuterium in the sample.

Solution:

In the ICON NMR configuration select the solvent 'None' in Lock/Shim Options → Solvent/Probe Dependencies.

Set Associated Shim Routine to XAU topshimfindsolvent2 (AU program from Bruker)

Set Associated Lock Routine to LOCK-OFF

What it does: The system does not try to lock to a deuterium signal. It then runs a ^1H spectrum, finds the biggest signal and uses that signal for 1D gradient shimming.

What it doesn't do: Optimise X-Y shims, lock (so there is no compensation for magnet drift or other magnetic field-related issues)

Tested on AVII+, AVIII running TS3.5, TS3.6