



A fresh approach in NMR spectrometry
for demanding applications

Company Introduction

QUAD Systems 2024



Background

Who are we ?

⌚ QUAD is a Swiss based private company bringing decades of collective experience in developing innovative NMR products.

What are our objectives ?

- ⌚ With the NMR market currently dominated by one major player, QUAD will provide alternative options for all NMR users.
- ⌚ Our goal is to rejuvenate innovation in NMR development for the whole community.

What do we offer ?

⌚ Full systems

- ⌚ Magnets & Shim systems
- ⌚ Consoles & Automations (sample changer, cooling units)
- ⌚ Probes

⌚ Replacements

- ⌚ Consoles + probes **ON YOUR** magnets + shim systems
- ⌚ Consoles + probes + shim systems (+ automation) **ON YOUR** magnet




Innovative products and performance

Competitive and economic price structure

Applicative and service team available within different countries (EU, US, CN)

QUAD NMR Full system



Collaborations with more than 30 NMR groups mainly in EU and North America.

QUAD NMR Probes

Routine probes:

- ④ 2CH probes for 300...600 MHz:
 - BBO/BBFO & BBI configuration
 - Z gradient: 6 Gauss/A*cm

R&D specific probes examples:

- ④ Ultra-high spin lock probes > 250 kHz on ^1H and ^{19}F
- ④ 3CH probe for Rh applications: H / P / ^{103}Rh -D 5 mm Z gradient
- ④ 3CH applications X / F / H 5 mm Z gradient
- ④ 3CH application X,F / P / H 5 mm Z gradient
- ④ 2CH 20 mm probes for e.g. rat heart applications $[^{13}\text{C}, ^{31}\text{P}, ^{23}\text{Na}] / ^1\text{H}$ 5 mm Z gradient



Piezo Actuators for
Tuning & matching






High precision
transmit / receive coil

We are open to build your R&D probes !

QUAD NMR Probes

Illuminated probes example:

②CH probes with light irradiation: glass fibre or LED arrays integrated into probe for UV, vis light irradiation



QUAD NMR Console

- ⌚ Dedicated to liquid state NMR
- ⌚ 2CH or 3CH, 300....600 MHz
- ⌚ ^2H or ^{19}F lock
- ⌚ Power Amplifier 1H = 100 W / X = 300 W / Lock = 10 W
- ⌚ Gradient amplifier 10 A
- ⌚ Separated preamp, near to magnet for highest signal/noise
- ⌚ 64CH shim power supply up to 22 bits resolution
- ⌚ Shim coil with max 47 shims



QUAD NMR Magnet

- ④ 400 and 600 MHz magnets with 54 mm bore
- ④ Highest field homogeneity and low drift rate < 4 Hz/h
- ④ Actively shielded for ease of siting and operation
- ④ Anti-vibration legs
- ④ Optimized for single story labs (< 3 m ceiling height)
- ④ High thermal efficiency cryostat for lowest possible Helium consumption: roughly 1 year hold time (LHe)



QUAD Cooling units

QUAD Peltier Cooler

- ⌚ Temperature range: +5 °C ... +25 °C
- ⌚ Energy efficient Cooling unit
- ⌚ No air compressor, no vibrations

QUAD Heat Exchanger

- ⌚ Temperature range: below -100 °C

Air-Jet Crystal cooler (SP Industries)

- ⌚ Temperature range: -90 °C ... +100 °C



QUAD Software suite







- ⌚ Intuitive Java coded software and subdivided into panels (no command lines)
- ⌚ Graphical instrument setting (wiring, probes management, shim...)
- ⌚ Role and permission management to regulate access to functionalities
- ⌚ Preparation, acquisition, process and export workflow
- ⌚ Tool for sequence development & shape designing
- ⌚ Integrated pulse sequences library
- ⌚ Pulse solver to compute the power of the pulses



- ⌚ Java coded software
- ⌚ Automation management
- ⌚ User interface is easy to pick up
- ⌚ Priority concept, Day/Night management

QUAD Software suite: SPINit



QUAD Software suite: SPINPlanner

The screenshot displays the SPINPlanner software interface. At the top right, there is a 'Login' button and a 'Logout' button. The main area is divided into several sections:

- Position:** A circular NMR tube rack with 100 positions, numbered 1 to 100. Position 1 is highlighted in red.
- Study/Experiment Settings:** A panel on the right showing details for the selected experiment, including 'Study: DemoStudy', 'Sample: Cyclosporin', 'Solvent: Benzene', 'Priority: NORMAL', 'Lock: Automatic lock', and 'Shim: Automatic'. It also lists the 'Series' of experiments.
- Series List:** A table listing the experiments to be performed:

Series	Name	Value
Proton (PulseAcquisition)	AcqTimeMode	Time
id: 256	Acquisition	2.7330718208s
nucleus: 1H		
Carbon (JMOD)		
id: 257		
nucleus: 13C, 1H		
id: 258		
nucleus: 1H		
HSQC (HSQC_edetgppsp.3_2D)		
id: 259		
nucleus: 1H, 13C		
HMBC (HMBC_acbigpl2ndqf_2D)		
id: 260		
nucleus: 1H, 13C		

- Schedule:** A table at the bottom showing the overall schedule of experiments:

Status	Tube	User	ID	Serie	Study	Sample	Temper...	Tune	Spin	Lock	Shim	Scan	Total	Queue	Scheduled...	Completion
●	1	Demo	256	Proton	DemoStudy	Cyclosporin		⊙	○	○	○	1min36s	9min16s		11:00:15	11:09:31
●	1	Demo	257	Carbon	DemoStudy	Cyclosporin		⊙	○	○	○					
●	1	Demo	258	COSY	DemoStudy	Cyclosporin		⊙	○	○	○					
●	1	Demo	259	HSQC	DemoStudy	Cyclosporin		⊙	○	○	○					
●	1	Demo	260	HMBC	DemoStudy	Cyclosporin		⊙	○	○	○					

Three callout boxes are overlaid on the interface:

- Login:** A box pointing to the 'Login' button at the top right.
- Select a position:** A box pointing to the circular NMR tube rack.
- Prepare experiments:** A box pointing to the 'Series' list.
- Control schedule:** A box pointing to the 'Schedule' table at the bottom.

QUAD Library

- ④ 650+ pulse sequences
- ④ 480+ applications
 - Pulse sequence associated with predefined acquisition parameters and process list
- ④ Protocols
 - Predefined list of applications
- ④ Shapes pulses
 - 170 Tx shapes
 - 9 gradient shapes
- ④ Predefined process list (1D, 2D, diffusion, JRES...)

You don't need to search too hard...

550+ pulse sequences loaded & tested with more to come...

sales@quadsystems.tech
www.quadsystems.tech



A fresh approach in NMR spectrometry
for demanding applications

**Thank you for your
attention !**

