

## Useful Keyboard Shortcuts

Esc – Cancel current mode  
a – Assignment  
b – Baseline correction  
c – Crosshair, Shift toggles 1 & 2 crosshairs (homonuclear only)  
f – Full spectrum  
i – Integration  
j – Multiplet analysis  
k – Peak picking: Manual threshold  
l – Reference peak  
m – Manual Zoom  
p – Pan  
P – Phase mode  
r – Reference graphically  
v – Restore (cuts etc.)  
w – Apodization  
x – Cut (1D)  
z – Toggle Zoom (in) mode: 1x/2x/3x z = horiz./vert./box zoom  
Z – Zoom out

### Ctrl/Cmd +

d – Export PDF  
i – Open directory  
k – Peak picking: Peak by Peak  
n – New document

### Ctrl/Cmd +

o, s, p, x, c, v, a:  
The obvious, like in any other program...